

Coupled LC-GC for the Analysis of Pesticide in Food

Ryoichi Sasano, Yutaka Nakanishi
Saika Technological Institute Foundation
75-2 Kuroda, Wakayama-city, Japan

Key Words: LC-GC, SPE, large volume injection, pesticide, food

1. Introduction

Analysis of pesticide residues in food is not so easy to do because various foods contain complex matrix. Recently several reports have been described that on-line LC-GC system is useful for the analysis of complex samples such as pesticide residues in food. The high sample capacity and wide range of separation performance of LC can be utilized in selective cleanup of sample and GC has high separation efficiency. However, coupling of LC to GC is not a trivial matter because the introduction of a large amount of LC fraction into a GC column requires the use of special techniques to separate the solvent from the sample. One of the subject for coupling LC with GC is a technique of eliminating water which is contained in LC fraction. Another subject is a technique of transferring analytes to GC.

The purpose of this study is coupling reversed-phase HPLC with GC for attaining automated analysis of pesticide residues in food. An interface equipped with a SPE cartridge is developed for coupling LC with GC. Finally, the performance of this LC-(SPE)-GC system has been evaluated by inspecting the data with analyzing chlorpyrifos in ginger.

2. Experiment

The system is shown in the following. With adding water to the LC fraction, the diluted fraction is loaded on the SPE cartridge. The analyte is adsorbed on the SPE. The cartridge was dried with nitrogen gas. The analyte is eluted from the the SPE cartridge with 20 μ L of hexane by syringe pump and the eluate is directly injected into a GC injector via a needle. The GC injector consists of a PTV injector equipped with a stomach shaped insert. The stomach shaped insert makes it possible to accept all of the injected eluate from the SPE cartridge.




Fig. Schematic flow diagram of LC-GC system.

LC-GC System Conditions

HPLC (MIDAS;Spark, Agilent 1100)		Interface Injector (LaviStoma; EMINET)	
Injection:	100 μ L, Sample loop	Insert:	Stomach Type Insert
Column:	3.0 mm i.d. \times 100 mm Inertsil ODS-3	Solvent Vent:	10 sec, Purge flow 150 mL/min
Solvents:	A: Acetonitrile/Water (50/50) B: Acetonitrile	Splitless:	3 min
	Flow rate 0.5 mL/min	Inj. Temp.:	70°C(3min)-120°C/min-220°C/min (3min) -50°C/min-260°C(10min)
Detector:	UV 210 nm	GC/MS (QP-5050A; Shimadzu)	
Interface SPE (LGI-S100)		Column:	0.25 mm i.d. \times 30 m, 0.25 μ m Inert Cap 5MS/Sil
SPE:	2 mm i.d. \times 10 mm C18	Oven:	70°C(3min)-20°C/min-280°C(4min)
Diluting:	Water 0.4 mL/min	Carr. gas:	He, 1 mL/min
Purge:	N ₂ gas, 0.5 min	MS:	SCAN;50-550 m/z
Elution:	Hexane, 20 μ L		

Sample preparation

20 g portions of homogenized ginger were extracted with 70 mL of acetonitrile. The extract solution was adjusted to 100 mL with water.

Then 1 mL of the extract was cleaned up through an C18 cartridge (50 mg) for avoiding HPLC column from deteriorating, and adjusted to 2 mL with 70% acetonitrile-water for LC-(SPE)-GC/MS analysis.

3. Result and Discussion

The performance of the system was investigated with respect to the LC-(SPE)-GC process by a ginger spiked with 0.1 μ g/g of chlorpyrifos. Chlorpyrifos elute around 8.2min. The fraction eluting from 7.8 to 8.6 min was transferred to the SPE interface. The excellent chromatography obtained can be observed in Fig. 2. The cleaning achieved with the LC process is highly satisfactory.




Fig.2 HPLC chromatogram of a ginger spiked with 0.1 μ g/g of chlorpyrifos (a) and a standard solution of it (b). Marked fraction transferred to the SPE cartridge.




Fig.3 GC/MS-SCAN chromatogram of the LC-(SPE)-GC/MS analysis of a ginger spiked with 0.1 μ g/g of chlorpyrifos.

4. Conclusion

Coupling of LC with GC has been accomplished by a new SPE interface and a PTV injector equipped with a stomach shaped insert. The LC-(SPE)-GC system provides very high efficiency and selectivity performance, and then allows automated analysis of pesticide residues in food.