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
1. Formic acid TBDMS
  2. Acetic acid TBDMS
  3. Propionic acid TBDMS
  4. Isobutyric acid (Methyl propionic acid) TBDMS
  5. Butyric acid TBDMS
  6. 2-Methyl butyric acid TBDMS
  7. Isovaleric acid (3-Methyl butyric acid) TBDMS
  8. Valeric acid TBDMS (Pentanoic acid)
  9. 2-Methylvaleric acid (2-Methylpentanoic acid) TBDMS
  10. 3-Methylvaleric acid (3-Methylpentanoic acid) TBDMS
  11. Isocaproic acid TBDMS (4-Methylvaleric acid, 4-Methylpentanoic acid, Isohexanoic acid) TBDMS
  12. Caproic acid (Hexanoic acid) TBDMS
  13. Lactic acid 2 TBDMS
  14. Succinic acid 2 TBDMS
1. 2-Hydroxybutanoic acid (2 TBDMS)
  2. 2-Hydroxy-3-methylbutyric acid (2 TBDMS)
  3. 2-Hydroxyisocaproic acid (2 TBDMS)
  4. Glyceric acid (3 TBDMS)
  5. 3-Phenylactic acid (2 TBDMS)
  6. 2-Hydroxyglutaric acid (3 TBDMS)
  7. Orotic acid (3 TBDMS)
  8. 4-Hydroxyphenylactic acid (3 TBDMS)
  9. Citric acid (4 TBDMS)

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化合物名	%RSD
Formic acid TBDMS*	1.7
Acetic acid TBDMS*	4.2
Propionic acid TBDMS	2.0
Isobutyric acid TBDMS	1.7
Butyric acid TBDMS	1.0
2-Methylbutanoic acid TBDMS	2.5
Isovaleric acid TBDMS	1.7
Valeric acid TBDMS	2.2
2-Methylvaleric acid TBDMS	2.1
3-Methylvaleric acid TBDMS	3.1
Isocaproic acid TBDMS	2.9
Caproic acid TBDMS	2.6
Lactic acid 2 TBDMS*	1.7
Succinic acid 2 TBDMS	1.1

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化合物名	%RSD	化合物名	%RSD
Formic acid TBDMS	0.8	2-Hydroxybutanoic acid 2 TBDMS	-
Acetic acid TBDMS	2.4	2-Hydroxy-3-methylbutyric acid 2 TBDMS	5.7
Propionic acid TBDMS	2.2	2-Hydroxyisocaproic acid 2 TBDMS	7.6
Isobutyric acid TBDMS	2.3	Glyceric acid 3 TBDMS	18.4
Butyric acid TBDMS	2.8	3-Phenylactic acid 2 TBDMS	12.8
2-Methylbutanoic acid TBDMS	1.4	2-Hydroxyglutaric acid 3 TBDMS	6.3
Isovaleric acid TBDMS	2.8	Palmitic acid TBDMS	8.1
Valeric acid TBDMS	2.5	Orotic acid 3 TBDMS	1
2-Methylvaleric acid TBDMS	1.5	Oleic acid TBDMS	4.1
3-Methylvaleric acid TBDMS	1.7	4-Hydroxyphenylactic acid TBDMS	8.9
Isocaproic acid TBDMS	1.2	Citric acid 4 TBDMS	2.9
Caproic acid	1.8		
Lactic acid 2 TBDMS	1.6		
Succinic acid 2 TBDMS	0.3		

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