

Development of an analytical method for volatiles by solid phase adsorption-analytical derivatization-solvent elution technique using online SPE-GC/MS system

OShoya Arakawa, Tomonori Asai, Ryoichi Sasano (AiSTI SCIENCE Co., Ltd.)

1. Solid phase adsorption-analytical derivatization-solvent elution technique






Fig 1. Overview of online SPE-GC/MS system LVI-S250 and SPL-M100 (AiSTI SCIENCE)

Headspace volatiles are collected and derivatized in a SPE cartridge, the entire eluate is injected into a large volume injection port.

2. Experimental method

Pretreatment workflow of volatile analysis (trimethylsilylation)




GC-MS condition

Methods	①Volatiles (non-derivatization)	②Volatiles (trimethylsilylation)
PTV Injector	LVI-S250 (AiSTI Science)	LVI-S250 (AiSTI Science)
Injector temp.	70°C(0.3 min)-100°C/min-290°C	150°C(0.5 min)-100°C/min-290°C
GC-MS	7890B + 5977B (Agilent Technologies)	7890B + 5977B (Agilent Technologies)
Split	100 mL/min (70 kPa, 0.3 min) → Splitless (3 min) → 50 mL/min	1:25
Flow mode	Constant flow, 1.4 mL/min	Constant flow, 1.0 mL/min
Pre-column	0.25 mm i.d. × 0.5 m	0.25 mm i.d. × 0.5 m
Column	DB-WAX, 0.25 mm i.d. × 60 m, df=0.25 μm	VF-5ms, 0.25 mm i.d. × 30 m, df=0.25 μm
Oven temp.	40°C (4 min) - 5°C/min - 230°C (10 min), Total 52.5 min	80°C (4 min) - 10°C/min - 210°C - 20°C/min - 310°C (2 min), Total 24 min
Transferline temp.	250°C	290°C
Ion source temp.	260°C	260°C
Acquisition mode	Scan (m/z 41-400)	Scan (m/z 41-500)

3-1. Repeatability and Linearity by musty odor


Repeatability (non-derivatization)

Sample: Spiked ultrapure water
Concentration: 10 μg/L
Sample volume: 5 mL
Additive: NaCl 2 g
Sampling volume: 20 mL
Flow : 150 μL/sec
Injection mode: Large volume injection (ⓐ)
GC Oven temp.: 40°C (4 min) - 20°C/min - 230°C (10 min)




Linearity (non-derivatization)

Sample: Spiked ultrapure water
Concentrations: 0.4, 2, 4, 6, 8 μg/L
Sample volume: 5 mL
Additive: NaCl 2 g
Drawing volume: 20 mL
Flow : 150 μL/sec
Injection mode: Large volume injection (ⓐ)
GC Oven temp.: 40°C (4 min) - 20°C/min - 230°C (10 min)




The system repeatability and sampling linearity were confirmed


3-2. Applicability and annotated volatiles



3-3. Discrimination test by coffee samples



In non-derivatization method, only MCB showed characteristic group of compounds



Derivatization method showed that pyrazines and phenols contributed

4. Discussions and Conclusion

Advantages of the derivatization method

- Analysis can be performed using the same configuration for metabolome analysis
- Shorter analysis time and high-throughput

Disadvantages of the derivatization method

- Detection of >C3 components is difficult (elution is faster than reagent)
- Sensitivity is lower compared to the non-derivatization method + large-volume injection

Future issues

- Establishment of a database of derivatized volatiles
- Simultaneous analysis with taste components (amino acids, organic acids, sugars, etc.)