

Simultaneous Multi-Component Analysis of PFAS in River Water Using an Online SPE-LC/MS System

Introduction

Perfluoroalkyl and polyfluoroalkyl substances (collectively referred to as "PFAS") are widely used in daily life as processing aids for fluoropolymers, in paints, water repellents, emulsifiers, firefighting agents, and nonstick cookware. However, concerns have been reported regarding environmental pollution and various toxic effects on humans due to their non-volatility and persistence, which cause them to remain in the environment for extended periods.

Currently, perfluorooctanoic acid (PFOA), perfluorooctane sulfonic acid (PFOS), and perfluorohexane sulfonic acid (PFHxS) are targeted for testing in water and soil. However, it is anticipated that the number of substances subject to regulation will increase in the future.


This application introduces a multi-component simultaneous analysis method for PFAS in river water using the SPL-W100 online SPE-LC/MS system, which performs fully automated solid-phase extraction and measurement simply by loading the sample.

Overview of the Online SPE-LC/MS System


This system achieves "reduced sample volume," "shortened pretreatment time," and "online operation" through the use of the dedicated solid-phase cartridge Flash-SPE and the mixed injection valve system [MiVS]. Furthermore, overlapping the pretreatment (solid-phase extraction) and instrument measurement enables more efficient analysis.



Online SPE-LC/MS System
SPL-W100



High-throughput analysis with overlapping preprocessing and measurement






Mixed Injection Valve System [MiVS] (Patent)

This valve features an optimized flow path design. It mixes the eluent from the solid phase with a diluent (such as water) within the valve, then switches the valve to introduce the mixture into the LC column. For example, in reverse-phase mode analysis, it can be used to mix with water to reduce the solvent ratio, adjust the pH of the eluent, or facilitate derivatization.

Dedicated Solid-Phase Cartridge Flash-SPE (Patent)

Since the filling amount is only a few milligrams, small-scale solid-phase extraction is achievable, and it is also effective for reducing solvent usage.



SPL-W100 for SPE-LC system

Sample



Information

This application is the result of a joint research project with the Saitama Prefectural Environmental Science International Center (Mr. Takamine and Mrs. Takazawa).

Key Word

PFAS
Solid-phase extraction
Online SPE-LC/
MS/MS

AiSTI SCIENCE

Product

SPL-W100
Flash-SPE WAXs





AiSTI SCIENCE CO.,Ltd.

TEL : +81-73-475-0033

E-Mail : as@aisti.co.jp

HP : www.aisti.co.jp

Experimental Methods



Measurement Conditions

[Device]

SPL-W100(AiSTI SCIENCE)、LCMS-8045(Shimadzu)

[LC Conditions]

Delay Column: Inertsil ODS-3, 3 μ m, 3.0 mm ID $\times$ 33 mm (UP)Analytical Column: Inertsil ODS-3, 3 μ m, 2.1 mm ID $\times$ 75 mm

Mobile Phase A: 2 mM ammonium acetate-water

B: 2 mM ammonium acetate-MeOH-acetonitrile (1/1)

Flow rate: 0.3 mL/min

Gradient: B.Conc. 20% (0 min) $\rightarrow$ 30% (1-2.5 min) $\rightarrow$ 40% (3-3.5 min) $\rightarrow$ 100% (11-15 min) $\rightarrow$ 20% (16 min)Column temperature: 40 $^{\circ}$ C

[MS Conditions]

Ionization mode : ESI Negative

Measurement Mode : MRM



Result

Addition recovery tests were conducted using tap water, mineral water, and river water. The recovery rate was calculated from the following equation, using samples added to purified water as the absolute calibration curve (1 point). For the three components PFOSA, NMePFOSA, and NePFOSA, recovery was not obtained. This is because they lack anionic charge, resulting in weak retention on the WAX, and because loading with hydrophobic organic solvents likely prevents retention on the WAX. For the other components, good recovery rates and reproducibility were obtained for all samples.

$$\text{Recovery Rate (\%)} = (\text{Added Area Value} - \text{Blank Area Value}) / (\text{Purified Water Added Area Value} - \text{Blank Area Value}) \times 100$$


Table: Recovery Rates and Reproducibility of PFAS in Each Sample

No.	成分名	Tap water		Mineral water		River water	
		20 ppt		20 ppt		20 ppt	
		Addition n=5	RSD (%)	Addition n=5	RSD (%)	Addition n=5	RSD (%)
1	PFBA	101	4.0	96	2.4	96	4.2
2	PFMPA	92	3.4	92	2.7	85	2.9
3	PFPeA	92	3.7	89	2.6	91	5.9
4	PFMBA	88	5.5	81	1.8	75	3.8
5	PFBS	99	3.3	103	4.5	99	6.0
6	4:2 FTS	99	5.3	94	4.9	84	7.9
7	PFHxA	96	3.5	96	2.7	97	3.4
8	PFEESA	97	3.0	99	4.0	96	3.7
9	HFPO-DA	99	6.6	103	2.1	103	4.2
10	PFPeS	99	6.0	101	2.8	102	4.2
11	PFHpA	97	2.6	95	2.2	97	0.9
12	PFHxS	99	3.9	103	4.2	103	3.1
13	6:2FTSA	96	3.4	97	4.5	99	2.7
14	PFOA	91	6.5	100	3.0	100	3.3
15	PFHpS	100	6.2	102	6.8	105	5.0
16	8:2 FTUCA	99	3.8	102	2.8	104	1.6
17	PFNA	90	3.4	93	2.7	94	2.6
18	PFOS	101	6.8	101	2.9	104	7.2
19	8:2FTSA	94	3.0	98	4.3	98	1.5
20	PFDA	90	3.9	97	3.2	93	2.4
21	FOSA	N.D.	—	N.D.	—	N.D.	—
22	NMeFOSAA	104	7.2	108	4.2	108	6.5
23	NeFOSAA	101	5.2	106	3.3	108	3.9
24	PFUnA	89	2.2	97	2.2	100	2.6
25	PFDS	101	7.4	102	2.4	109	5.5
26	PFDoA	86	1.9	103	2.7	104	5.3
27	NMeFOSA	N.D.	—	N.D.	—	N.D.	—
28	NeFOSA	N.D.	—	N.D.	—	N.D.	—
29	PFTTrDA	92	2.7	116	1.6	117	3.0
30	PFTeDA	92	5.3	118	1.0	121	4.3
31	PFHxDA	83	2.8	103	1.5	106	7.3
32	8:2 diPAP	96	4.3	113	2.6	118	2.4
33	PFOcDA	80	6.7	91	5.2	89	3.1


※Calculations are based on absolute area values without applying

[Unrecovered components]

FOSA



NMeFOSA



NeFOSA

