# Development of a Rapid, Simple, and Highly Sensitive Multi-Component Simultaneous Analysis Method Combining Liquid-Liquid Extraction and Solid-Phase Extraction for PFAS in Food and Serum

Ryoichi Sasano<sup>1</sup>, Heesoo Eun<sup>2</sup> (1 AiSTI SCIENCE Co., Ltd., 2 National Agriculture and Food Research Organization (NARO))

✉ E-mail: sasano@aisti.co.jp

## Summary

In the analysis of PFAS in food and serum, a large number of matrix components complicate sample preparation, requiring significant time and effort. Moreover, as the number of analytes increases, pretreatment becomes even more challenging. In this study, highly polar ionic matrix components were first removed by liquid-liquid extraction (LLE) to reduce the load on the solid phase. Subsequently, PFAS were retained and concentrated using anion-exchange solid-phase extraction (SPE), which simultaneously eliminated nonionic matrix components to improve both sensitivity and purification. Furthermore, an online SPE-LC/MS system was employed to automate the process from solid-phase extraction to LC/MS measurement, enabling rapid, simple, and highly sensitive multi-component analysis.

## Experimental Methods

A PFAS mixed standard solution containing 35 compounds (FUJIFILM Wako Pure Chemical Corporation) was used. Samples included mackerel (fish), milk, and human serum. The pretreatment flow is shown in Figure 1. Each sample was extracted and deproteinized with acetonitrile (ACN), followed by liquid-liquid partitioning with a hydrophobic solvent. The organic extract was then introduced into an online SPE-LC/MS system composed of an SPL-W100 online SPE unit (AiSTI Science) coupled with an LCMS-8045 triple quadrupole mass spectrometer (Shimadzu Corporation). A WAXs cartridge (Flash-SPE WAXs, 3 mg) packed with anion-exchange resin was used as the solid-phase column.


Figure 1. Pretreatment flow for PFAS analysis using liquid-liquid extraction and online solid-phase extraction.



Figure 2. Pretreatment procedure for PFAS analysis


Online SPE-LC/MS System

### Ion-Exchange Solid-Phase Extraction.

Anionic PFAS in the ACN-TBME organic phase were retained on the positively charged WAXs solid-phase cartridge through ion-exchange interaction.

PFAS were eluted by alkaline 0.1% ammonia in methanol-water, which neutralized the WAXs functional groups and released the analytes.

The process of online SPE-LC/MS

## Results and Discussion

### Deproteinization

Since food and serum samples contain abundant proteins that can interfere with solid-phase performance, deproteinization was carried out in the initial extraction. Acetonitrile was selected as the solvent after comparing its compatibility with subsequent liquid-liquid partitioning. The acetonitrile-to-water ratio was set at 2:1.

### Liquid-Liquid Extraction

After deproteinization with acetonitrile, TBME and HCl were added to the supernatant to perform LLE. This transferred PFAS into the organic phase while removing polar matrix components in the aqueous layer. The addition of HCl promoted transfer of polar PFAS by reducing ionization.

### Solid-Phase Extraction

The hydrophobic organic extract was loaded onto the WAXs anion-exchange cartridge. Ionic PFAS were retained, while nonionic interferences were removed during washing. Although retention of target analytes can be affected by coexisting anions such as nitrate, preliminary removal of these species via LLE maintained strong retention. Increasing the sample loading volume further improved sensitivity.

### Online Automation

By introducing the entire eluate from the solid phase directly into the LC/MS system, sensitivity was significantly improved. This enhancement allowed for the use of smaller sample volumes while maintaining reliable quantification.

### Analysis Time

The entire process from solid-phase extraction to LC injection was fully automated, requiring only 12 minutes per operation. The LC/MS analytical cycle time was 15 minutes, enabling continuous and parallel operation to obtain measurement results every 15 minutes.

### Blank Control

To minimize contamination, all solvents used in the solid-phase extraction process were pre-filtered through an adsorbent material integrated into the pretreatment system. This setup effectively reduced background signals and ensured high analytical reliability.

### Recovery Test


Recovery and reproducibility (n = 5) were evaluated by comparing absolute peak areas of samples spiked with PFAS standard mixtures, unspiked samples, and standard solutions measured with the online SPE-LC system (Table 1).

Hydrophobic and nonionic compounds such as FOSA, NMeFOSA, and NETFOSA were not recovered, likely because they were not retained on the solid phase and eluted with organic solvents during washing. Reproducibility was satisfactory, with relative standard deviations (RSD) below 10% for most analytes.

## Conclusion

By combining deproteinization, liquid-liquid partitioning, and anion-exchange solid-phase extraction, this study achieved rapid, simple, and highly sensitive multi-component PFAS analysis even for complex matrices such as food and serum.


### Liquid-Liquid Extraction Method


For the partitioning solvent, tert-butyl methyl ether (TBME) and ethyl acetate—both immiscible with water but miscible with ACN—were evaluated. TBME showed the highest recovery efficiency and was selected as the optimal solvent.

### Effect of HCl Addition

Addition of HCl improved the recovery of hydrophilic PFAS such as PFBA and PFMPA by suppressing carboxyl group dissociation, increasing hydrophobicity and facilitating transfer to the organic phase.


### Evaluation of Liquid-Liquid Extraction

In serum samples, comparison between samples with and without LLE revealed that when the aqueous phase (ACN-water ratio 1:9) was directly loaded onto the WAXs column, recoveries were very low, suggesting that matrix interference inhibited ion exchange. In contrast, a 4:1 ACN-water ratio allowed PFAS retention through hydrophobic interactions, although recovery of highly hydrophobic PFAS slightly decreased, which could be improved by optimizing solvent ratios.

### Absolute calibration curve using the online SPE-LC/MS system.

### Recovery and Reproducibility Test

No.	Compound	Mackerel 2 ppb		Milk 1 ppb		Human Serum 1 ppb	
		Rec. %	RSD n=5, %	Rec. %	RSD n=5, %	Rec. %	RSD n=5, %
1	PFBA	71	7.1	82	3.4	106	5.1
2	PFMPA	69	5.9	88	4.1	91	1.8
3	PFPeA	90	3.2	89	6.3	102	9.4
4	PFMBA	73	8.4	88	3.8	91	6.3
5	PFBS	107	2.4	79	6.5	108	11.1
6	4:2 FTS	94	8.1	60	6.0	121	8.0
7	NFDHA	147	4.3	94	4.9	185	7.3
8	PFHxA	93	3.5	84	4.2	94	6.4
9	PFEESA	90	4.0	89	2.1	105	5.7
10	HFPO-DA	87	4.5	85	6.3	91	8.8
11	PFPeS	107	4.1	84	7.2	106	9.2
12	PFHpA	94	3.8	85	3.6	95	5.7
13	PFHxA	106	4.3	92	8.6	107	6.1
14	6:2FTSA	130	4.7	81	6.3	192	5.1
15	PFOA	108	2.1	85	3.0	102	3.9
16	PFHpS	107	4.3	88	2.9	102	10.7
17	8:2 FTUCA	123	1.6	86	4.2	115	7.1
18	PFNA	101	1.6	88	4.3	102	4.7
19	FOSA	31	26.0	43	5.2	49	7.0
20	PFOS	96	3.7	80	3.8	82	5.1
21	8:2FTSA	142	8.1	80	8.1	133	10.7
22	PFDA	99	4.4	88	2.9	106	6.6
23	NMeFOSAA	113	5.1	97	4.1	135	8.2
24	PFNS	80	5.3	88	3.6	96	9.7
25	N-MeFOSA	-	-	-	-	-	-
26	NETFOSAA	51	4.9	87	7.5	93	7.4
27	PFUnA	49	3.2	85	2.6	84	4.7
28	PFDS	93	5.6	93	4.4	98	2.6
29	N-EtFOSA	-	-	-	-	-	-
30	PFDoDA	63	3.9	83	3.6	77	5.2
31	PFTrDA	74	2.6	80	5.0	87	7.5
32	PFTeDA	89	4.2	80	3.7	98	5.8
33	PFHxDA	74	6.0	86	3.3	98	5.8
34	8:2 diPAP	89	2.0	87	6.7	79	8.9
35	PFOcDA	113	2.6	94	3.4	138	5.3

✉ The recovery rate was calculated from the absolute peak area values (without using stable isotopes).

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