

# Evaluation of a fully automated analytical method for PFAS in environmental water

Mari Takazawa<sup>1</sup>, Ryoichi Sasano<sup>\*2</sup>, Shusuke Takemine<sup>1</sup>, Mamoru Motegi<sup>1</sup>, Koji Kosaka<sup>3</sup>

<sup>1</sup>Center for Environmental Science in Saitama, <sup>2</sup>AIsti SCIENCE Co., Ltd., <sup>3</sup>National Institute for Environmental Studies

E-mail: takazawa.mari@pref.saitama.lg.jp

## Introduction

Perfluoroalkyl substances (PFAS) are regarded as a persistent environmental concern on a global scale. In recent years, the technology of fully automated solid-phase extraction (online SPE) has attracted attention due to its potential to reduce costs, save time, and minimize environmental impact. In a previous study, a direct injection was attempted into HLB or WAX (Takazawa, et.al, 2024). However, the method proved inadequate for application to the sample containing high-ion contaminants. The present study proposes a methodology for the removal of ionic contaminants from a sample prior to its loading onto WAX.

## Material and Method

The automated analytical method is shown in Fig. 1. Depending on the analytes, several methods were used. These included up to 47 analytes, such as 13 PFCAs, 10 PFSA, and 24 other PFASs. The analytes were covered by ISO 21675 and EPA methods 533, 537, and 1633. Five picograms of the standard mixture was added to 1 mL of ultrapure water and used as test water. In the case of surface water test, PFOA and PFHxS were not spiked because their concentrations in the surface water were sufficient concentrations for this evaluation. In the 47 PFAS recovery tests, 10 pg of the standard mixture was added; however, the lower limit of detection has not been examined. Concurrent surrogate recovery tests were conducted; however, no surrogate correction was performed to properly evaluate the chemical behavior of each compound.



Online SPE-LC-MS/MS system





Fig. 1 Extraction procedures for PFAS with Online SPE-LC-MS/MS system.

## Results and Discussion

### PFOS · PFOA · PFHxS


#### Hydrophobic interaction solid-phase extraction.

##### Load and Wash



PFAS in water are retained on the HLB via hydrophobic interactions. Whereas hydrophilic matrix components are not retained and are removed.

##### Elution:




The retained PFAS are eluted from the HLB using methanol and collected. Whereas hydrophobic matrix components remained in the cartridge.

### 47 PFAS

#### Ion-Exchange Solid-Phase Extraction.

##### Load and Wash




Anionic PFAS in the organic phase were retained on the positively charged WAXs through ion-exchange interaction. Whereas non-anionic matrix components are not retained and are removed.

##### Elution:




PFAS were eluted by alkaline ammonia in methanol, which neutralized the WAXs functional groups and released the analytes.

#### Effect of Liquid-Liquid Extraction




For an organic phase, tert-butyl methyl ether (TBME) and ethyl acetate were selected, which were completely separated with a sample water. Ethyl acetate-TBME(1/1) showed the highest recovery efficiency and was selected as the optimal solvent.

#### The conventional method



#### The present method



#### Effect of HCl Addition

Addition of HCl improved the recovery of hydrophilic PFAS such as PFBA and PFMPA by suppressing carboxyl group dissociation, increasing hydrophobicity and facilitating transfer to the organic phase.





Fig. 6 Recoveries of PFAS in surface water with SPE Cartridges of WAX.




Fig. 2 Pretreatment procedure for PFAS analysis.





Fig. 3 The process of online SPE-LC-MS/MS.

#### Dilution effect in the elution process via online SPE-LC-MS/MS

##### No dilution



##### Dilution with water




Fig. 4 Peak shape improvement by sample dilution with water after extraction with online SPE system.





Fig. 5 Absolute calibration curve using the online SPE-LC-MS/MS system.

#### Table 2 Recovery and Reproducibility Test

No.	Analyte	Operation blank pg/mL	River water pg/mL	River water + native standard RSD, n=6	
				%	Rec.
1	PFBA	0.2	4.1	3.2	92
2	PFPeA	0.2	2.0	1.8	95
3	PFHxA	1.0	3.6	4.2	96
4	PFHpA	0.2	1.4	3.4	99
5	PFOA	0.3	4.2	2.9	100
6	PFNA	0.2	1.2	3.0	103
7	PFDA	0.2	0.2	4.3	102
8	PFUnA	0.2	0.2	2.8	98
9	PFDoDA	0.1	0.1	3.6	95
10	PFTriDA	0.2	0.0	2.6	97
11	PFTeDA	0.1	0.0	3.3	96
12	PFHxDA	0.1	0.0	5.1	89
13	PFOcDA	0.2	0.0	12.6	91
14	8:2 diPAP	N.D.	N.D.	3.2	89
15	PFBS	N.D.	1.4	4.3	93
16	PFPeS	N.D.	N.D.	3.0	94
17	PFHxS	N.D.	0.2	9.3	87
18	PFHpS	N.D.	N.D.	3.2	89
19	PFOS	N.D.	0.5	7.1	84
20	PFNS	N.D.	N.D.	4.0	97
21	PFDS	N.D.	N.D.	8.0	109
22	8:2 FTUCA	0.1	0.0	2.4	110
23	4:2 FTS	N.D.	N.D.	18.7	102
24	6:2FTSA	N.D.	N.D.	5.6	104
25	8:2FTSA	N.D.	N.D.	7.8	103
26	PFMPA	N.D.	N.D.	2.4	98
27	PFMBA	N.D.	N.D.	4.5	97
28	NFDHA	N.D.	N.D.	3.0	88
29	HFPO-DA	N.D.	1.8	5.0	104
30	PFEESA	N.D.	N.D.	2.4	99
31	NMeFOSAA	N.D.	N.D.	4.0	104
32	NEtFOSAA	N.D.	N.D.	3.3	110

Sample: Brackish water (Wakayama prefecture, Japan)  
\* The recovery rate was calculated from the absolute peak area values (without using stable isotopes).

#### The volume of extract loaded onto the solid phase



#### Optimization of elution volume from SPE cartridge

